

Small Drug Safe – MSD-1 (Tiny)

Schedule 8 Drug Safe

The first Schedule 8 compliant Drug Safe specifically aimed at low volume storage for use by mobile vets and medicine-dependent households.

Features:

- Fully Welded Steel Construction
- Steel thickness 10mm body 10mm door
- Lock: 6 lever pick-resistant key lock with 2 keys
- Anti-tamper mechanisms:
 - Internal rear protection lock plate
 - Welded hinge protector

Configuration:

- No shelves
- 4 rear fixing holes – vertically orientated to facilitate fixing to studs
- 4 base fixing holes – for floor fixing

Fixing: It is a requirement of compliance that an S8 drug safe is fixed to either a wall or floor through the holes provided in back and rear of the safe. The surface being fixed to will determine the type of fixing (coach screws, dynabolts, anchor bolts, etc. These are not included)

Maintenance:

The Tiny requires very little maintenance except preventing the build-up of dust which may penetrate the locking mechanism over time. To clean the safe, just wipe with a non-abrasive cloth using a soft household cleaner such as 'Spray and Wipe.'



****Now with ergonomic handle for opening and closing**



Model	Width (mm)	Depth (mm)	Height (mm)	Weight (Kg)
External	190	190	270	20
Internal	170	170	250	